
MONARK
EXERCISE AB

MANUAL
English/Svensk

Ergomedic 891 E

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Congratulations on your new exercise device.

Designed by MONARK EXERCISE AB, Sweden.

Monark has been the world's leading manufacturer of high ergometers and exercise cycles for more than 40 years.

Contents	Page
Assembly instructions	3 - 4
Battery replacement	7
Brake belt – brake surface	8
Brake belt – replacement/adjustment	7 - 8
Chain	9
Crank bearing	11
Digital meter	6
Flywheel bearing	10
General	3
General about exercise	11 - 12
Operation instruction	5 - 7
Service checklist	10 - 11
Warrenty	3

GENERAL

It is important that you keep your exercise device clean and properly lubricated. Most important is to protect the chromed and zined parts, but also painted parts benefit from the same protection. When cleaning and lubricating, be sure to check that all screws and nuts are properly tightened. Be sure that all moving parts as crank and flywheel is working normal and that no unnormal play or sound exists. I.e. play in bearings causes fast wearing and with that follows a highly reduced lifetime.

PLEASE NOTE: The production number of your exercise device is placed according to fig 3, page 4.

WARRENTY

As on any quality product there may be an exceptional fault due to material or manufacture. If such a fault should arise on your exercise cycle, please return to the place of purchase for necessary repair. Monark products and parts are guaranteed against defects in materials and workmanship for a period of one year from the initial date of purchase of the unit.

Parts found to need replacement due to normal wear and tear, such as brake belts, are not covered. This guarantee covers parts only, not labor costs associated with the repair.

This guarantee does not apply to cases of abuse or vandalism, nor does it extend to any injury or loss to person or property caused directly or indirectly by any Monark products.


In the event of a defect in material or workmanship during the warranty period, Monark Exercise will repair or replace (at its option) the product. Monark Exercise will do so at its expense for the cost of materials but not for labour or shipping

PARTS BELOW ARE NOT ASSEMBLED

Assemble crank handle marked R (Right) on to the right hand side. The crank handle axel has a right hand thread and must be threaded onto the crank clockwise.

Assemble crank handle marked L (Left) on to the left hand side. The crank handle axel has a left hand thread and must be threaded onto the crank counter clockwise.


Untighten the two nuts on the frame for equipment. mount the weight holder and remount the nuts tightly.

ASSEMBLY


Assemble the stand

Assemble first the support stand as shown to the right.
Mount all the screws and nuts
before you tightend them.

Mount the device on the stand

Assemble first the device to the
help frame with four locking nuts
and two screws (MVBF M8) as
shown to the right.
Mount the help frame to the
support stand with four screws,
washer and nuts.

OPERATING INSTRUCTIONS

NOTE!

When using the Cardio Rehab 891E severe physical tension can occur. Therefore it is recommended that persons, who are not used to fitness training or don't feel well, contact a doctor/physician for consultations.

Monark Exercise Ergometer device Model 891E is a test device, applied with a braking system, where the brake force can be adjusted with weights on a weight carrier.

Weights are available in 0,5 kg and 0,1 kg.

Note: 0,1 kg is the minimum weight, which is the weight of the weight carrier. The weight carrier can also be placed in an upper level and then will no braking force be applied. Disengage by pushing the knob. See fig below.

The Cardio Rehab, model 891E, is equipped with an electronic meter showing pedal revolutions per minute (RPM), exercise time in minutes and seconds (TIME), total pedal revolutions (COUNT). CAL is not a valid figure.

When pedaling the subject stores energy in the flywheel. This is braked by means of a brake belt which runs around the bigger part of the brake of the flywheel. The workload is changed either by using other pedaling speed or increasing or decreasing the tension of the brake belt against the flywheel by means of the workload tension device.

NOTE!

At transport the tension device should be somewhat tightened to prevent the brake belt from falling off the flywheel.


Fig. 4.

OPERATING INSTRUCTIONS

Digital meter

Features

Time	Timer	up/down
RPM	Pedal	rev./min
COUNT	Pedal	rev. total
CAL	Approx. Cal. rel. to RPM	


Fig. 5.

891-007

Keys and functions

1. MODE key

Use the MODE-key to cycle through the functions. To activate the automatic SCAN function use MODE to advance to SCAN function. SCAN will appear on the display.

When scan function has been activated, each function will be displayed in 5 seconds intervals throughout your exercise session.

2. SET key

Pressing this key will flash to show computer is in setting mode for programming your EXERCISE TIME. If you hold down this key for 2 seconds you can advance the function value at a faster rate.

Press Mode key to advance to TIME functions, and use SET to enter your desired time.

Each press of SET will advance time by one minute. The time set symbol ((o)) will appear on the display. When pedaling is started, measurement of remaining exercise time will begin and the unit will count down to 0:00 at which time beeper will sound for 5 seconds.

3. RESET key

The RESET key will clear the values for TIME, COUNT and CAL individually.

NOTE:

1. Display will return to normal when not pressed any key for 5 seconds, or trainer is not used.
2. The computer starts automatically when one of the keys is pressed or when the trainer is used, that is meter gets rpm indication.
3. Auto shut off function when not pressed any key or device not used for 4 minutes.

CALORIE function


Since the meter is not connected to any force sensor, the calorie function is just intended to be used as a guide at different rotation speed.

OPERATING INSTRUCTIONS

Replacement of batteries

Pull the connector out of the housing and take the lid panel of.

The batteries, 2 x 1,5 V, size AA (R6), which are placed in a holder on the backsize of the panel, can easily be changed. See fig. 6.

Replacement of brake belt/cord

Removed the front cover over the weight carrier by losen the screws on each side of the cover.

Take the return strap, see fig. 8, 9, and lift up the weight basket until it locks in its upper postion.

Loosen the cord/belt bracket, see fig. 9, and take away the brake from the tension hub.


Loosen or cut away the knot or tie up the knot at the other end of the belt. After that take the belt from the device.

When assembling a new brake cord/belt, first enter one end into the hole in the tension hub, see fig. 7, and tie a knot and let the knot fall into the bigger part of the hole.

Then assemble the new belt exactly as the old one.

NOTE: When replacing the brake cord/belt it is recommended to clean the brake surface.

See BRAKE SURFACE - BRAKE BELT.

Adjusting the brake cord/belt tension

At first check that the brake belt is lying correct on the flywheel surface.


Put 1 kg in the weightt basket. Rotate the flywheel by hand. the basket shall now lift up so the distance to the flywheel is at least 40 mm, and maximum 80 mm. If this is the case, the brake belt has to be loosened or tightened a little at the tension center hub. If the basket i too low shorten the belt somewhat and if the basket is too high leghten the cord somewhat.

Lock the weight basket in its upper postion and after that loosen the cord bracket somewhat so that the cord length can be adjusted. Tighten the bracket again, release and check if the measurment above ar OK when the flywheel is rotated by hand.

Repeat the above if necessary.

NOTE: The basket will always give the correct workload wherever it hangs between the upper stop position and the lower stop just above the flywheel. The measurments above just gives a good margin - in each direction

Brake belt contact surface - brake belt

The brake belt should be checked now and then to ensure that it has not suffered excessive wear. If it looks worn it should be replaced.

Deposits of dirt on the brake belt and on the contact surface may cause the unit to operate unevenly and will also wear out the brake belt. The brake belt contact of the flywheel surface should then be ground off with a fine sandpaper and any dust removed with a clean dry cloth.

Dismantle cover, see "Replacement brake belt".

Set the weight basket in its upper position.


Take off the brake belt to the side.

Grind with a fine sandpaper, see fig. 10.

Grinding is easier to perform if a second individual cautiously and carefully pedals the device.

Irregularities on the brake belt contact surface are removed by means of a fine sandpaper or an abrasive cloth. Otherwise unnecessary wear on the brake occurs and the unit can become noisy.

Always keep the brake belt contact surface clean and dry. No lubricant is allowed to be used.

Chain 1/2" x 1/8"

It is strongly recommended that a chain solvent be used to keep the chain clean. Excess dirt built up on the chain will cause excess wear. A chain lubricant and solvent for normal road bikes may be used.

Check the lubrication and tension of the chain at regular intervals. In the middle of its free length the chain should have a minimum play of 5 mm. See fig 11. When the play in the chain is about 20 mm (<1 inch) the chain must be tightened otherwise it will cause unnormal wear of the chain and chainwheels. Because of this it is always recommended to keep the chain play as little as possible. When the chain has become so long that it can no longer be tightened with the chain adjusters it is worn out and shall be replaced with a new one.

ADJUSTING CHAIN

Remove left and right frame cover. See "REPLACEMENT of BRAKE BELT".

To adjust the chain the hub nuts should be loosened. Loosening or tightening the nuts on the chain adjusters will then move the hub and axle forward or backward. Adjust according to above recommendation. Then tighten the nuts on the hub axle again. See fig 13.



CHAIN REPLACEMENT

Loosen the chain adjusters as much as possible. Dismantle the chain lock and remove the chain. Put on a new chain and assemble the chain lock. The spring of the chain lock should be assembled with the closed end in the movement direction of the chain. Use a pair of tongs for dismantling and assembling the spring. See fig 12.

Adjust chain adjusters to chain play according to above. Tighten axle nuts firmly.

Put on frame covers again.

NOTE: At assembly the flywheel has to be parallel with the center line of the frame. Otherwise the chain and chainwheels makes a lot of noise and wears out very rapidly.

Service and maintenance

REPLACEMENT OF THE FREEWHEELING SPROCKET

Remove left and right frame cover. See "REPLACEMENT OF BRAKE BELT".

Dismantle the chain as described on page 9.



Loosen the axle nuts and lift off the flywheel. Remove the axle nut, washer, chain adjuster and spacer on the freewheel side. Place the special remover (part No. 9100-14) in the adapter and place the spacer and axle nut outside. See fig 14.

NOTE: Do not tighten the axle nut completely. It must be possible to loosen the adapter-sprocket half a turn.

Replace sprocket-adapter and assemble the new parts in reverse order according to the above.

LUBRICATION SPROCKET

The sprocket should be lubricated with a few drops of oil once a year. Incline the cycle somewhat to make it easier for the oil to reach the bearing. See fig 15.

Service checklist

CRANK BEARING

The crank bearing is long term greased and needs normally no supplementary lubrication. If problem arises, please contact your Monark dealer.

FLYWHEEL BEARING

The bearings in the flywheel are lifetime greased and require normally no maintenance. *If problem arises, please contact your Monark dealer.*

SERVICE CHECKLIST

Check the:

- ☐ chain is snug and there is no play on the pedal crank
- ☐ pedal crank is secure to the crank axle
- ☐ pedals moving smoothly, and is the pedal axle clear of dirt and fibres
- ☐ pedals are securely fitted to the pedal crank
- ☐ flywheel rotating smoothly and central
- ☐ handlebars screws lubricated
- ☐ brake belt does not show significant signs of wear
- ☐ pedals and chain are lubricated

THE IMPORTANCE OF REGULAR EXERCISE

The human body is built for action – not for rest. Once upon a time this was a necessity: the struggle for survival demanded good physical condition. But optimal function can only be achieved by regularly exposing the heart, circulation, muscles, tendons, skeleton and nervous system to some loading, i.e. training.

In the old days the body got its exercise both in work and at leisure. In our modern society, however, machines have taken over an ever increasing share of the tasks which were formerly accomplished with muscular power alone. Our life has at an accelerated tempo been dominated by sitting, riding and lying. Thus, the natural and vital stimulation that tissues and internal organs receive through physical exercise has largely disappeared. Certain tissues such as muscles, bone and blood and also a number of bodily functions can adapt to inactivity – and to stress. Studies have proved that if you use 30 minutes for exercise like brisk walking, running, bicycling, swimming or skiing 2-3 times a week, your condition has been improved by some 15 per cent after a few months. The efficiency of the heart muscle will increase and joints and muscles grow in strength. The capillary density increases in the trained muscle and their enzymatic activities are enhanced. The body adapts to the new demands. The perceived exertion at a given rate of exercise becomes reduced.

With increased physical activity fatness is concentrated, the appetite functions "safer", you can eat more without risk for overweight and thereby the risk of lack of important essential food nutrients decreases. For many individuals the effect of habitual physical activity also improves the wellbeing and it is a good feeling to have a potential to cope with straining situations.

What kind of exercise to choose?

1. You should have fun when exercising. Choose something you find pleasure in doing regularly.
2. To get a good effect out of the training you should choose a form of exercise that engages large muscle groups. Then the demand of increased blood flow and oxygen transport will be so great that heart will increase its pump capacity. Jogging, calisthenics, aerobic dancing, bicycling, swimming, skiing and walking are excellent examples of exercises meeting this requirement.

IN A FEW MONTHS YOU CAN GET 10-15 YEARS YOUNGER

If you cycle 30 minutes a few times a week you can lower your condition age with 10-15 years! Scientifically this is described as a reduction on the biological age. Externally, you are your usual self. Internally, however, you feel much younger. In other words: You can work harder. You feel more alert and healthy. Your ability to handle stress and problems increases. There are few better ways to improve your physical condition than to cycle. It does not over-tax your joints. It builds up your condition progressively and at your own pace – and you can make your training fit weather conditions.

DO I LOOSE WEIGHT WHEN I CYCLING?

Yes! You do lose calories. A few miles on your bike every day over one year, you will have lost the equivalent of 20 pounds of body fat. You will achieve best results if you combine exercise with healthier eating. A little less sugar, less butter on your bread or less fat in your frying pan. And a few miles on your bike every day. In a year you will have lost 20 pounds.

DO I GET STRONGER?

Cycling strengthens the muscles of the back, abdomen and legs. Daily chores become easier. Cycling also makes your heart stronger. Your pulse rate gets lower even when you exert yourself a little extra. Regular exercise also has a favourable influence on high blood pressures.

HOW DO I TRAIN?

1. Warm up 3-5 minutes with a low pedal resistance. Pedal about 12 mph (20 km/h).
2. Increase the resistance until you feel the training "somewhat hard". Keep the speed for 2-5 minutes. Get off the Ergometer and rest a few minutes. Cycle again and then rest. Train at your own pace and with a comfortable pedal resistance. After a few weeks you can increase the resistance.
3. Before ending, pedal a few minutes with a light resistance, in order to step down your training.

Total time about 30 minutes.

Strength training:

1. Give yourself a thorough warm-up.
2. Pedal with a heavy resistance for 5-10 seconds, then rest 45-60 seconds. Repeat this 5-10 times.

It is a good idea to combine your cycle training with gymnastics for 5 minutes, as this will give you a physiologically well-balanced form of training.

(Elderly people and physically weak persons should consult a doctor before starting their training.)

Notes:

Lined area for notes, consisting of multiple horizontal dotted lines.

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