IMPORTANT NOTE: PLEASE READ THE FOLLOWING PRODUCT INFORMATION CAREFULLY BEFORE USING THE LACTATE SCOUT ANALYZER!

Product information about the Lactate Scout handheld analyzer with Lactate Scout test strips to determine lactate in whole blood


SET MODE (entered by switching on the analyzer by scroll button / without inserted test strip)

The Lactate Scout can be switched on either by inserting a test strip (see „Display 1“) or by pressing the scroll button for 1 sec. If you switch on by using the scroll button, you will access the „set mode“ which allows you to enter the test strip code, to check the memory, to launch the function control, timer functions and some other options of the analyzer. **Turning the scroll wheel upwards or downwards, the available options will be flashed in succession on the display. Press briefly on the scroll button when you have reached the option you want to access.** Depending on the chosen option, you can look for stored data, change current settings, or perform a device and sensor check. If an option is selected by pressing the scroll button, you can change its existing settings by turning and pressing the scroll wheel, e.g. sound volume, duration length, code characters etc. **By pressing the button briefly, the displayed setting will be confirmed and stored.**

After confirmation of an individual setting, the next setting related to the previous appears automatically (so after setting the hours, you have to set the minutes etc.) until the available settings of this option are completed. **After completing all positions of an option, you will be back in the general setting mode again,** to select another option, to switch off the analyzer or to switch in the measurement mode by inserting a test strip.

Pressing the scroll button 3 sec. or more, „OFF“ appears on the display and the analyzer turns off. If the device is in set mode and is not being used for more than 3 minutes, it turns off automatically. If the Lactate Scout is turned off while an option has been chosen but not confirmed/stored yet, the settings or changes may not be saved. This also happens when a test strip is inserted at this situation: The analyzer changes immediately to the measurement mode then (see „Display 1“). Please always ensure that you have confirmed/stored your settings before.

set code: Presets the device to the two characters of the given test strip calibration code (0 1 2 3 4 5 6 7 8 9 R b. c. d. E. F as shown on the label of the test strip vial) in two steps.

mem: To review the stored readings in the internal memory. The last result will be displayed at first, added by the corresponding memory ID (001 - 250). After 250 measurements, the first readings will be replaced by the last (e.g. measurement no. 256 will be stored at pos. 006). After 1 sec., the memory ID disappears again and the display shows the related date and time of this measurement below its displayed result. Measurement readings taken in the step test mode are shown with a flashing set of squares to indicate the step test level (PRL flashes on left side, Ld middle, AFL right).

dSc: To check the functionality of the Lactate Scout test strips and analyzer with the Lactate Scout test solution. Because the test solution is different in its characteristics to human blood, the accuracy of its readings will be lower (values can differ from 8.9 to 11.1 mmol/l) than readings of whole blood. **No calibration required! Control or calibration solutions from 3rd parties (e.g. for lab devices) will not work properly with the Lactate Scout system.**

Selecting dSc option, „ctrl d+s“ will be displayed fixed, and „dSc“ (device + sensor check) will flash. Insert a new test strip, the code appears for 1 sec (ensure correct code setting!) and a flashing drop appears. For 20 sec., signalled by beeps, the analyzer is waiting to adsorb a sample of control solution as usual. Open the bottle of test solution, and press a drop out of a clean surface (e.g. finger nail). Lower the analyzer with the inserted test strip into the drop then. After sample uptake or after the 20 sec., the function control starts like a regular measurement. After 10 seconds, a control value appears briefly, followed by „OK“ if no problem was detected. Wipe the remaining drop away. If something was wrong, „Er 4“ will display. Repeat it with a new test strip or contact customer service.

Snd: Sound volume is shown by the length of the bar underneath the letters „Snd“ (full or half sound volume). If no bar is shown, the analyzer is on mute: Please note, you get no acoustic signals (beeps) if then.

dAt: Setting of day, month, hour, minute and year (please note this sequel!).

dur: Setting the duration of an exercise/load phase in minutes and seconds. These figures also function as the starting time for the countdown („cdn“). The maximum duration time is 10 minutes.

clC: This function works like a simple stop watch, measuring the time in seconds after pressing the scroll button.

cdn: Countdown (related to „dur“ function), helps to control the time of exercise/load phase to the next measurement.

If you feel unfamiliar in using the scroll wheel button: Take care that the scroll button is no longer turned when pressing, what could cause the selection of another option than chosen originally. By continual pressing, the analyzer switches off.

POWER SUPPLY / BATTERIES

To operate the Lactate Scout analyzer, use two standard batteries type LR03/AAA (1.5 V). **Never plug an external power supply or 3rd party components into the socket of the analyzer! This socket is only for connection with the Lactate Scout PC products. Disregarding this limitation may occur several damages, and the warranty will be void!**

The operational capability using alkaline batteries as delivered with is at least one year. Carbon-Zinc batteries can be used also, but do not mix two different kinds of batteries at the same time! Additionally to the „bAt“ message, there is a battery symbol which appears on the top left of the display as a warning that the batteries will have to be changed soon. At this point, there is still enough power to run the Lactate Scout analyzer for some days. The data storage and settings will not be lost during the replacement of the batteries.

The Lactate Scout analyzer must be switched off when replacing the batteries. Gently push up the battery cover on the bottom. Take out old batteries and replace with new ones. **Please make sure that you follow the polarity markings to orient the batteries properly (+ / -)!** Push the lid back with care until it is firmly closed. The Lactate Scout analyzer is now ready to use again.

Do not throw used batteries into the domestic waste and ask your local dealer for proper dispose!

Distributed by:

EKF-diagnostic sales GmbH
Ebendorfer Chaussee 3
D-39179 Barleben/Magdeburg
Tel.: +49 (0)39 203-785 0
Fax: +49 (0)39 203-785 24
www.ekf-diagnostic.com
sales@ekf-diagnostic.com

LIFE SUPPORT POLICY


LACTATE SCOUT - PRODUCTS ARE NOT AUTHORIZED FOR USE AS CRITICAL COMPONENTS IN LIFE SUPPORT DEVICES OR SYSTEMS WITHOUT THE EXPRESSLY WRITTEN APPROVAL OF THE MANUFACTURERS PRESIDENT OF SENSLAB GMBH, GERMANY.

FDA owner / operator ID number: 064575
FDA Device Listing Number: E 184489
The device is 510(k) exempt under 21CFR 862.1450

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IMPORTANT NOTE: PLEASE READ THE FOLLOWING PRODUCT INFORMATION CAREFULLY BEFORE USING THE LACTATE SCOUT ANALYZER!

Lactate **SCOUT** - analyzer

Product information about the Lactate Scout handheld analyzer with Lactate Scout test strips to determine lactate in whole blood

INTENDED USE

The Lactate Scout analyzer is to be used only with the Lactate Scout test strips to determine the lactate content in capillary blood. For in-vitro-application only! The Lactate Scout analyzer and the supplied test strips can be used both by individuals and by health personnel.

WARNING NOTICES

Please ensure the Lactate Scout analyzer and accessories are kept strictly away from children, to prevent accidents caused by improper handling and to avoid the risk of small parts being swallowed! If you get results which are unexpectedly high or low, the measurement may be affected by substances like sweat or alcohol from the disinfecting swab. If sweat or alcohol is mixed with the blood sample, the reading may be incorrect. Clean the area with soap and water and dry thoroughly. If you think the measurement is wrong, repeat the measurement with a new test strip and check the meter with the Lactate Scout test solution as described under the menu display "dSC". **Read the instruction manual and package insert before using the Lactate Scout!**

Measuring lactate carries the risk associated with exposure to blood and infection: Please make sure items containing blood such as test strips, lancets and alcohol swabs are disposed of properly according to local regulations, to avoid any risk to third parties.

Always wear gloves when you are carrying out a measurement on another person!

CLEANING

The device has to be cleaned by wet tissue, soaked with dish detergents and water. Keep liquids from entering the slot for test strips, PC port and scroll button. **If the Lactate Scout analyzer, the lancet device or any other items came in contact with blood, use suitable tools as alcohol swabs or disinfection sprays and wear disposable gloves.** To avoid damages of the display of the Lactate Scout analyzer, dry and clean it immediately after disinfection. All other items should dry at room temperature.


Medical personnel should observe the internal regulations of their particular institution.

STORAGE AND HANDLING

- Like any other electronic devices, the Lactate Scout analyzer should be handled and used carefully.
- Protect from moisture, permanent sun light and temperatures lower than 0°C and higher than +50°C.
- Rough handling such as dropping or pushing the analyzer can lead to restricted ability of operation. If happened, contact customer service before using it again!
- To avoid battery leakage, remove them if the device will not be used for a considerable time.
- In cases of mistakes in operating, error messages, or failure see notes in the instruction manual.


APPLICATION / USAGE

- **The analyzer has only to be used with the Lactate Scout test strips, follow their product insert.**
- Make sure that the test strip you insert comes from a vial with the same code that appears on the display of the analyzer. If the code does not match, reset the code number of the analyzer as described in the instruction manual.
- ~~Clean the fingertip or earlobe carefully with water and disinfect the area with alcohol swab.~~
- Remove any residue from the skin with a clean tissue and dry it thoroughly before lancing/puncturing.
- The surface of the lancing device must be clean and disinfected. Use a new lancet for each test!
- Puncture the fingertip or earlobe with the lancing device.
- Squeeze slightly the surrounding area to get the blood sample.
- To avoid contamination, wipe away the first drop of blood. Use the next drop that forms for measurement.
- Lower the tip of the inserted test strip to the drop of blood until it is absorbed, confirmed by a signal.
- After the test, take the test strip in the middle, remove it from analyzer and dispose of it properly according the local regulations. Do the same with the used lancet and any other disposable materials that have contacted blood.

DISPLAY 1 (after insertion of the test strip, measuring mode)

- Device switches on, the test strip code on the analyzer appears for 2 sec.
- If the code does not match that of the vial of the test strip inserted, pull out the test strip, and the analyzer will switch off. For code setting switch on the analyzer again and follow the instructions below.
- If the code flashes after inserting the test strip then the analyzer has recognized that 24 measurements were done using test strips with an identical code. To avoid wrong measuring values the analyzer asks for confirmation or correction of the code. For confirmation press briefly on the scroll button for each of the two digits.
- The display will be replaced by a flashing drop symbol, indicating readiness.
- By turning the scroll button, options of the step test mode can be selected ("PrL": pre-load/resting value, "Ld": main/exercise load, "AFL": after load/recovery value).
- No display: After 3 min. without any activities, the display switches off to save battery power. Press the scroll button to activate again.

DISPLAY 2 (after sample uptake by test strip)

- When the required sample volume is taken up completely, it will be confirmed by a beep signal.
- The display shows "LAC", and a bar line shows the time left (10 sec. measurement duration).
- If readings are above 25 mmol/l "Hi"(high)" appears, and "Lo" if lower than 0.5 mmol/l.
- In the case of these error messages pull out the test strip and use a new one.
- After successful measurement the result will be displayed, signalled by another beep.
- The results are automatically stored with time, date, measurement mode and memory ID.
- To delete/reject results (only possible during the test strip is inserted continually) turn the scroll button to "del." (delete) and press briefly. The display shows "- -" as confirmation.
- After 3 min. without any activities after inserting the test strip, the display switches off.
- When the test strip is removed, the device switches off completely.


DISPLAY 3 (switching on by scroll button / without inserted test strip: Set mode)

- The analyzer can be switched on by holding the scroll button down for 1 second.
- If you press the scroll button longer than 3 sec., the device turns off.
- Switching on by scroll button, "set code" is the first option that will appear.
- Turning the scroll wheel will show other options: To select an option, press the scroll button briefly:

set code: Setting of test strip code	Snd: Sound volume	cl: Stop watch
mem: Displays stored readings	dAt: Date and time	cdn: Countdown/Timer
dSC: Device and sensor control	dur: Duration of step length (Details given in the manual)	

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Lactate SCOUT
 (UK)
 Product information about the Lactate Scout test strips for one-time disposable use with the Lactate Scout and Lactate Scout+ analyzer to determine lactate in whole blood.
 Please read carefully before using!

INTENDED USE
 Lactate Scout test strips are only to be used with the Lactate Scout analyzer to determine the lactate concentration in capillary blood. For in-vitro-application only! The Lactate Scout test strips and the Lactate Scout analyzer can be used both by individuals and by health personnel.

WARNING NOTICES
 Read the instruction manual and package insert before using the Lactate Scout! Avoid any risk of test strips being swallowed! Keep away from children! Always wear gloves when you are carrying out a measurement on another person! The handling with blood carries generally the risk of infections. Please make sure items containing or coated with blood (such as test strips, lancets and alcohol swabs) are disposed of properly according to the local regulations, to avoid potential abuse, contamination and other risks to third parties.

STORAGE
 Keep the test strips protected from moisture, light, and temperatures above +25°C/77°F! Keep refrigerated! Transport and store the test strips only in their original vial! Take out the required test strips for immediate use only and close the vial again, to keep the other test strips protected continually. Before using refrigerated test strips, take the vial out and keep it at regular temperature for 20 min at least. When used, analyzer and test strips should be at same temperature. After first opening the vial of test strips, the storage stability is three months or until expiry date if it is sooner. Afterwards, inaccurate measurement results will occur increasingly.

WRONG READINGS/INTERFERENCES
 Wrong little rests of sweat, saliva, desinfection or alcohol on the skin surface may cause problems of sample uptake or leads to wrong readings (as too high results)! Ensure that the area (fingertip or earlobe) is cleaned carefully by water and dried before each sampling again. Hurry up with the sample uptake after puncturing: Only fresh capillary blood can be used for meaningful results!
 Also taken drugs (as paracetamol) as well as illness, infection, diet, fat meals, physical or psychical stress may influence the measurement results differently.

HANDLING (ONE-TIME DISPOSABLE USE)
 After cleaning and puncturing of the sampling area, a new, non-used Lactate Scout test strip has to be inserted into the Lactate Scout analyzer. Therefore, hold the test strip in the middle and do not touch the lip! The hands should be clean and dry. By inserting the test strip, the calibration code will be displayed for 1 s. Ensure that it is identical to the code as printed on the currently used vial of test strips. Also control the expiry date!

Lower the Lactate Scout analyzer with the tip of the inserted test strip from top down to the fresh drop of blood (if possible, wipe away the first drop before and take the second volume). Avoid smearing the blood sample, and do not touch the skin! The required sample volume (0.5 µl) will be uptake automatically. When the display shows LAC, the measurement runs and the result will be shown after 15/10 s (Lactate SCOUT / SCOUT+). By removing the test strip, the result will be stored automatically (incl. date, time and other data), and the Lactate Scout analyzer switches off.

Error messages: By readings above 25.0 mmol/l, "HI" appears, and "Lo" by readings lower than 0.5 mmol/l. As explained in the manual, if the temperature is lower than +5°C/41°F or higher than +45°C/113°F, a warning symbol appears and no measurement occurs.
Health note: There are no known health risk from using Lactate Scout test strips according to the instructions. SensLab products are not authorized for use as critical component in life support devices or systems! Persons prone to excessive bleeding or infection, or persons feeling ill have to consult a doctor generally before puncturing and sampling.

TECHNICAL DATA
 Measurement method: Enzymatic-amprometric detection of lactate
 Reagent per test strip: Lactate oxidase, electron mediator, additives
 Coefficient of variation: 3.0 - 8.0% (depends on concentr.; validated average)
 Sample volume: 0.5 µl
 Measuring range: 0.5 - 25.0 mmol/l

TEST SOLUTION (option):
 Contains buffer solution, additives and liquid with a defined lactate concentration (see imprint). The test solution is suitable only for a general function control (without calibration)! To use within 90 days after first opening and before expiry date. Other calibration solutions, liquids or control materials cannot be used properly with a Lactate Scout!

For additional information visit: www.lactatescout.com

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